

Medical body balanced care tape

192


Feature

- Strong adhesive force
- Using soft materials
- Through humidity by special adhesive materials.

Application

Medical body balanced care tape

Structure



Properties

general properties

Item	Adhesive force (N/inch)
192	19

※Measure method Adhesive force JIS Z0237
PET#25 backing, peeling speed 300mm/min, 180° angle.

Precautions on use

All technical data are prepared based on the tests and measured values carried out in the laboratory of KGK Chemical Corp. as the standard.
However, product characteristics may vary greatly depending on environment and adherend.
Therefore, regarding these characteristic data, it is a reference value, not a guaranteed value.
Before using it please make sure that this product is suitable for the intended use and environment.

Caution on storage

Please be sure to put it in a box and keep it.
Please choose a cold and dark place not to be exposed to direct sunlight for the storage location.
In particular, please do not expose to high temperature and high humidity

(temperature 30 °C or more and humidity 50% or more forbidden). 940 Minaminagai Tokorozawa-City saitama-Pref

The warranty period of the tape is unopened and it is six months after shipment.

KGK Chemical Corp.

359-0011 Japan

Tel : +81 4 2944 5151

Mail : postbox@kgk-tape.co.jp

URL : <http://www.kgk-tape.co.jp/>